

## Product datasheet for **HP208858**

### RNF41 Human qPCR Primer Pair (NM\_005785)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10193
Forward Sequence:	ATCGCAGAGCTGGAGAAGACGT
Reverse Sequence:	TGTCTCCTCCAGTTCTGAAGG
ACCN:	<a href="#">NM_005785</a> , <a href="#">NM_005785.2</a> , <a href="#">NM_005785.3</a> , <a href="#">BC032637</a> , <a href="#">BC032637.1</a> , <a href="#">BC011626</a> , <a href="#">BC024284</a> , <a href="#">BC036267</a> , <a href="#">BC095418</a> , <a href="#">BI562407</a> , <a href="#">BP201682</a> , <a href="#">BU170070</a> , <a href="#">NM_005785.4</a>
UniProt ID:	<a href="#">Q9H4P4</a>
Synonyms:	FLRF; NRDP1; SBBI03
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)