

Product datasheet for **HP208870**

USPL1 Human qPCR Primer Pair (NM_005800)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10208
Forward Sequence:	AATGTCATCCCTGAGTGGCACC
Reverse Sequence:	GCAAGTCATTCTGTGTAAGCCT
ACCN:	<u>NM_005800</u> , <u>NM_005800.1</u> , <u>NM_005800.2</u> , <u>NM_005800.3</u> , <u>NM_005800.4</u> , <u>BC038103</u> , <u>BC038103.2</u> , <u>BC012594</u> , <u>BC036027</u> , <u>BI752017</u>
UniProt ID:	<u>Q5W0Q7</u>
Synonyms:	bA121O19.1; C13orf22; D13S106E
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)