

## Product datasheet for **HP208891**

### AP3S2 Human qPCR Primer Pair (NM\_005829)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10239
Forward Sequence:	GATTGAGGCTCAAAACAGGCTGG
Reverse Sequence:	CAAACGGGACAGGTTGGGAAC
ACCN:	<a href="#">NM_005829</a> , <a href="#">NM_005829.1</a> , <a href="#">NM_005829.2</a> , <a href="#">NM_005829.3</a> , <a href="#">NM_005829.4</a> , <a href="#">BC002785</a> , <a href="#">BC002785.2</a> , <a href="#">BC007773</a> , <a href="#">BC010020</a> , <a href="#">BC010615</a> , <a href="#">BX098014</a> , <a href="#">NM_005829.5</a>
UniProt ID:	<a href="#">P59780</a>
Synonyms:	AP3S3; sigma3b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)