

Product datasheet for **HP208895**

p40 (RABEPK) Human qPCR Primer Pair (NM_005833)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10244
Forward Sequence:	CACATCATCGGCAGCCATTGGA
Reverse Sequence:	GGAGGATTCCAAGTGTCTCTGG
ACCN:	<u>BC000503</u> , <u>NM_005833</u> , <u>NM_005833.1</u> , <u>NM_005833.2</u> , <u>NM_005833.3</u> , <u>BC000503.2</u> , <u>BC053541</u> , <u>BC053541.1</u> , <u>BC047023</u> , <u>BC065725</u> , <u>NM_005833.4</u>
UniProt ID:	<u>Q7Z6M1</u>
Synonyms:	bA65N13.1; p40; RAB9P40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)