

## Product datasheet for **HP208911**

### CDK2AP2 Human qPCR Primer Pair (NM\_005851)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10263
Forward Sequence:	CGCTCCTTTCAGACCGCTGTTT
Reverse Sequence:	CCATCTCCTCTATGACTGACAGC
ACCN:	<a href="#">NM_005851</a> , <a href="#">NM_005851.1</a> , <a href="#">NM_005851.2</a> , <a href="#">NM_005851.3</a> , <a href="#">NM_005851.4</a> , <a href="#">BC002850</a> , <a href="#">BC002850.2</a> , <a href="#">BC016704</a> , <a href="#">BI763791</a> , <a href="#">BM559572</a> , <a href="#">BQ016039</a> , <a href="#">BT006909</a> , <a href="#">NM_005851.5</a>
UniProt ID:	<a href="#">O75956</a>
Synonyms:	DOC-1R; p14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)