

Product datasheet for **HP208937**

MARCHF6 Human qPCR Primer Pair (NM_005885)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10299
Forward Sequence:	GCTCACTACTGACGCTGCCATT
Reverse Sequence:	CTCTCAACCACACCAGGCTGAT
ACCN:	<u>NM_005885</u> , <u>NM_005885.1</u> , <u>NM_005885.2</u> , <u>NM_005885.3</u> , <u>BC013761</u> , <u>BC039293</u> , <u>BC046148</u> , <u>BC111494</u> , <u>BC136461</u> , <u>BC136462</u> , <u>BC142679</u> , <u>BC142694</u> , <u>BX116291</u> , <u>NM_005885.4</u>
UniProt ID:	<u>O60337</u>
Synonyms:	DOA10; FAME3; MARCH-VI; MARCH6; RNF176; TEB4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)