

Product datasheet for **HP208959**

MCM4 Human qPCR Primer Pair (NM_005914)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4173
Forward Sequence:	CTTGCTTCAGCCTTGGCTCCAA
Reverse Sequence:	GTCGCCACACAGCAAGATGTTG
ACCN:	NM_005914 , NM_005914.1 , NM_005914.2 , NM_005914.3 , BC031061 , BC031061.1 , BF793151 , BM781972 , BQ058022 , BQ689703 , NM_005914.4
UniProt ID:	P33991
Synonyms:	CDC21; CDC54; hCdc21; IMD54; NKCD; NKGCD; P1-CDC21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)