

Product datasheet for **HP208983**

MMP16 Human qPCR Primer Pair (NM_005941)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4325
Forward Sequence:	GATTCAGCCATTTGGTGGGAGG
Reverse Sequence:	CCCTTCCAGACTGTGATTGGC
ACCN:	<u>NM_005941</u> , <u>NM_005941.1</u> , <u>NM_005941.2</u> , <u>NM_005941.3</u> , <u>NM_005941.4</u> , <u>BC069500</u> , <u>BC069500.1</u> , <u>BC075004</u> , <u>BC075005</u> , <u>NM_005941.5</u>
UniProt ID:	<u>P51512</u>
Synonyms:	C8orf57; MMP-X2; MT-MMP2; MT-MMP3; MT3-MMP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)