

## Product datasheet for **HP208993**

### MTHFD1 Human qPCR Primer Pair (NM\_005956)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4522
Forward Sequence:	TTGGACAGGCTCCAACGGAGAA
Reverse Sequence:	AGAAGTGGTGAGAGCCAGGACA
ACCN:	<u><a href="#">NM_005956</a></u> , <u><a href="#">NM_005956.1</a></u> , <u><a href="#">NM_005956.2</a></u> , <u><a href="#">NM_005956.3</a></u> , <u><a href="#">BC009806</a></u> , <u><a href="#">BC009806.2</a></u> , <u><a href="#">BC001014</a></u> , <u><a href="#">BC050420</a></u> , <u><a href="#">BQ959308</a></u>
UniProt ID:	<u><a href="#">P11586</a></u>
Synonyms:	CIMAH; MTHFC; MTHFD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)