

Product datasheet for **HP208999**

MYH10 Human qPCR Primer Pair (NM_005964)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4628
Forward Sequence:	GCTGATGGCAACTCTCCGAAAC
Reverse Sequence:	CTTCCAGGACACCATTACAGCG
ACCN:	NM_005964 , NM_005964.1 , NM_005964.2 , NM_005964.3 , BC000280 , BC008968 , BC117690 , BC117691 , BC144668 , BC150634 , BE799557 , BG421687 , BG723653 , BU173337
UniProt ID:	P35580
Synonyms:	NMMHC-IIB; NMMHCB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)