

Product datasheet for **HP209007**

Papillary renal cell carcinoma (PRCC) Human qPCR Primer Pair (NM_005973)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5546
Forward Sequence:	AGATGGCAGCAGGTTCAAGTGG
Reverse Sequence:	CTGTGCAGGATAGTACCACCA
ACCN:	NM_005973 , NM_005973.1 , NM_005973.2 , NM_005973.3 , NM_005973.4 , BC010450 , BC010450.1 , BC004913 , BG746454 , BM423590 , BU942864
UniProt ID:	Q92733
Synonyms:	RCCP1; TPRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)