

Product datasheet for **HP209028**

alpha Tubulin (TUBA4A) Human qPCR Primer Pair (NM_006000)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7277
Forward Sequence:	GGCAAGGAGATCATTGACCCAG
Reverse Sequence:	CATCAGGAGTGAGGTGAAGCCA
ACCN:	<u>NM_006000</u> , <u>NM_006000.1</u> , <u>NM_006000.2</u> , <u>BC009238</u> , <u>BC009238.2</u> , <u>BC018854</u> , <u>BP346611</u> , <u>BT006731</u> , <u>NM_006000.3</u>
UniProt ID:	<u>P68366</u>
Synonyms:	ALS22; H2-ALPHA; TUBA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)