

## Product datasheet for **HP209032**

### UQCRH Human qPCR Primer Pair (NM\_006004)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7388
Forward Sequence:	GCTCTGTGATGAGCGTGTATCC
Reverse Sequence:	GTTGTTAAAGAGTTTGTGGCCAC
ACCN:	<a href="#">BC001426</a> , <a href="#">NM_006004</a> , <a href="#">NM_006004.1</a> , <a href="#">NM_006004.2</a> , <a href="#">NM_006004.3</a> , <a href="#">BC001426.2</a> , <a href="#">BC001934</a> , <a href="#">BC001934.1</a> , <a href="#">BC015177</a> , <a href="#">BC093060</a> , <a href="#">BC107703</a> , <a href="#">BG472687</a> , <a href="#">BP366403</a> , <a href="#">BT007070</a> , <a href="#">NM_006004.4</a>
UniProt ID:	<a href="#">P07919</a>
Synonyms:	QCR6; UQCR8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)