

Product datasheet for **HP209035**

ARMET (MANF) Human qPCR Primer Pair (NM_006010)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7873
Forward Sequence:	CAGCCACCAAAATCATCAATGAGG
Reverse Sequence:	TCCACTGTGCTCAGGTCGATCT
ACCN:	NM_006010 , NM_006010.1 , NM_006010.2 , NM_006010.3 , NM_006010.5 , BC007282 , BC033055 , BC113588 , BC113590 , BE514267 , BT007110 , NM_006010.6
UniProt ID:	P55145
Synonyms:	ARMET; ARP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)