

Product datasheet for **HP209039**

LAGE3 Human qPCR Primer Pair (NM_006014)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8270
Forward Sequence:	CAAAGGGTGGTTGGGAAGGATC
Reverse Sequence:	AGTTGATGACGGAAATTCGGAGC
ACCN:	<u>NM_006014</u> , <u>NM_006014.1</u> , <u>NM_006014.2</u> , <u>NM_006014.3</u> , <u>NM_006014.4</u> , <u>BC015744</u> , <u>BC015744.1</u> , <u>BC019012</u> , <u>BC062330</u> , <u>BQ649305</u>
UniProt ID:	<u>Q14657</u>
Synonyms:	CVG5; DXS9879E; DXS9951E; ESO3; GAMOS2; ITBA2; Pcc1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)