

## Product datasheet for **HP209043**

### HM74 (HCAR3) Human qPCR Primer Pair (NM\_006018)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8843
Forward Sequence:	CTGTTTCCACCTCAAGTCCTGG
Reverse Sequence:	CAGTCTGAACGCCGCACATAGT
ACCN:	<u><a href="#">NM_006018</a></u> , <u><a href="#">NM_006018.1</a></u> , <u><a href="#">NM_006018.2</a></u> , <u><a href="#">BC047891</a></u> , <u><a href="#">BC047891.1</a></u> , <u><a href="#">BC038955</a></u> , <u><a href="#">BI908162</a></u> , <u><a href="#">NM_006018.3</a></u>
UniProt ID:	<u><a href="#">P49019</a></u>
Synonyms:	GPR109B; HCA3; HM74; Puma-g; PUMAG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)