

## Product datasheet for **HP209044**

### TCIRG1 Human qPCR Primer Pair (NM\_006019)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10312
Forward Sequence:	CGGCTTCATCTACAACGAGTGC
Reverse Sequence:	TCCAGGTAAGCATCGTGTGCT
ACCN:	<u><a href="#">NM_006019</a></u> , <u><a href="#">NM_006019.1</a></u> , <u><a href="#">NM_006019.2</a></u> , <u><a href="#">NM_006019.3</a></u> , <u><a href="#">BC018133</a></u> , <u><a href="#">BC018133.1</a></u> , <u><a href="#">BC032465</a></u> , <u><a href="#">BC032465.1</a></u> , <u><a href="#">BX388184</a></u>
UniProt ID:	<u><a href="#">Q13488</a></u>
Synonyms:	a3; Atp6i; ATP6N1C; ATP6V0A3; OC-116kDa; OC116; OPTB1; Stv1; TIRC7; Vph1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)