

## Product datasheet for **HP209063**

### HDAC6 Human qPCR Primer Pair (NM\_006044)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10013
Forward Sequence:	GCCTCAATCACTGAGACCATCC
Reverse Sequence:	GGTGCCTTCTGGTGACCAACT
ACCN:	<a href="#">BC011498</a> , <a href="#">NM_006044</a> , <a href="#">NM_006044.1</a> , <a href="#">NM_006044.2</a> , <a href="#">NM_006044.3</a> , <a href="#">BC011498.1</a> , <a href="#">BC069243</a> , <a href="#">BC005872</a> , <a href="#">BC013219</a> , <a href="#">BC013737</a> , <a href="#">BT006649</a> , <a href="#">NM_006044.4</a>
UniProt ID:	<a href="#">Q9UBN7</a>
Synonyms:	CPBHM; HD6; JM21; PPP1R90
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)