

## Product datasheet for **HP209068**

### VPS26C Human qPCR Primer Pair (NM\_006052)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10311
Forward Sequence:	CCCATCTACATGGTCTTCCCTAG
Reverse Sequence:	GGAAGTTCTCCGTGATGAGGTG
ACCN:	<u><a href="#">NM_006052</a></u> , <u><a href="#">NM_006052.1</a></u> , <u><a href="#">BC110655</a></u> , <u><a href="#">BX648191</a></u> , <u><a href="#">NM_006052.2</a></u>
UniProt ID:	<u><a href="#">O14972</a></u>
Synonyms:	DCRA; DSCR3; DSCRA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)