

Product datasheet for **HP209104**

NDRG1 Human qPCR Primer Pair (NM_006096)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 10397 |
| Forward Sequence: | ATCACCCAGCACTTTGCCGTCT |
| Reverse Sequence: | GACTCCAGGAAGCATTTCAGCC |
| ACCN: | <u>NM_006096</u> , <u>NM_006096.1</u> , <u>NM_006096.2</u> , <u>NM_006096.3</u> , <u>BC003175</u> , <u>BC003175.1</u> , <u>BC006260</u> , <u>BX648361</u> , <u>NM_006096.4</u> |
| UniProt ID: | <u>Q92597</u> |
| Synonyms: | CAP43; CMT4D; DRG-1; DRG1; GC4; HMSNL; NDR1; NMSL; PROXY1; RIT42; RTP; TARG1; TDD5 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)