

## Product datasheet for **HP209117**

### Cyclophilin E (PPIE) Human qPCR Primer Pair (NM\_006112)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10450
Forward Sequence:	GGCGGTGATTTCAAAACCAC
Reverse Sequence:	CAGGAAGAACTGAGAGCCATTGG
ACCN:	<a href="#">NM_006112</a> , <a href="#">NM_006112.1</a> , <a href="#">NM_006112.2</a> , <a href="#">NM_006112.3</a> , <a href="#">BC008451</a> , <a href="#">BC004898</a> , <a href="#">BC107736</a> , <a href="#">BU160972</a> , <a href="#">BU179717</a> , <a href="#">BU196539</a> , <a href="#">BU689929</a> , <a href="#">NM_006112.4</a>
UniProt ID:	<a href="#">Q9UNP9</a>
Synonyms:	CYP-33; CYP33
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)