

## Product datasheet for **HP209132**

### GM CSF Receptor alpha (CSF2RA) Human qPCR Primer Pair (NM\_006140)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1438
Forward Sequence:	CCTGTCAGGATTAACGTCTCGC
Reverse Sequence:	CATTGCTGGGAGGTTGAATCG
ACCN:	<u><a href="#">NM_006140</a></u> , <u><a href="#">NM_006140.1</a></u> , <u><a href="#">NM_006140.2</a></u> , <u><a href="#">NM_006140.3</a></u> , <u><a href="#">NM_006140.4</a></u> , <u><a href="#">BC002635</a></u> , <u><a href="#">BC002635.1</a></u> , <u><a href="#">BC071835</a></u> , <u><a href="#">BI868829</a></u> , <u><a href="#">BM149332</a></u> , <u><a href="#">NM_006140.6</a></u>
UniProt ID:	<u><a href="#">P15509</a></u>
Synonyms:	alphaGMR; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMCSFR-alpha; GMR; GMR-alpha; SMDP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)