

Product datasheet for **HP209137**

Hsp40 (DNAJB1) Human qPCR Primer Pair (NM_006145)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3337
Forward Sequence:	AGTTCAAGGAGATCGCTGAGGC
Reverse Sequence:	GCTGAAAGAGGTACCATTGGCAC
ACCN:	<u>NM_006145</u> , <u>NM_006145.1</u> , <u>NM_006145.2</u> , <u>BC002352</u> , <u>BC002352.2</u> , <u>BC013157</u> , <u>BC019827</u> , <u>BM795601</u> , <u>NM_006145.3</u>
UniProt ID:	<u>P25685</u>
Synonyms:	Hdj1; Hsp40; HSPF1; RSPH16B; Sis1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)