

## Product datasheet for **HP209146**

### NEDD8 Human qPCR Primer Pair (NM\_006156)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4738
Forward Sequence:	CAGAGGCTCATCTACAGTGGCA
Reverse Sequence:	GTCCATCACTGCCTAAGACCAC
ACCN:	<u><a href="#">NM_006156</a></u> , <u><a href="#">NM_006156.1</a></u> , <u><a href="#">BC104664</a></u> , <u><a href="#">BC038290</a></u> , <u><a href="#">BC064899</a></u> , <u><a href="#">BC104200</a></u> , <u><a href="#">BC104201</a></u> , <u><a href="#">BE548741</a></u> , <u><a href="#">NM_006156.3</a></u>
UniProt ID:	<u><a href="#">Q15843</a></u>
Synonyms:	NEDD-8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online >](#)