

## Product datasheet for **HP209154**

### **Nrf2 (NFE2L2) Human qPCR Primer Pair (NM\_006164)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	4780
<b>Forward Sequence:</b>	CACATCCAGTCAGAAACCCAGTGG
<b>Reverse Sequence:</b>	GGAATGTCTGCGCCAAAAGCTG
<b>ACCN:</b>	<a href="#">NM_006164</a> , <a href="#">NM_006164.1</a> , <a href="#">NM_006164.2</a> , <a href="#">NM_006164.3</a> , <a href="#">NM_006164.4</a> , <a href="#">BC011558</a> , <a href="#">BC011558.1</a> , <a href="#">BI910808</a> , <a href="#">BP228042</a> , <a href="#">BX649047</a> , <a href="#">BX649071</a> , <a href="#">NM_006164.5</a>
<b>UniProt ID:</b>	<a href="#">Q16236</a>
<b>Synonyms:</b>	HEBP1; IMDDHH; Nrf-2; NRF2
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)