

Product datasheet for **HP209171**

Nurr1 (NR4A2) Human qPCR Primer Pair (NM_006186)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4929
Forward Sequence:	AAACTGCCCAGTGGACAAGCGT
Reverse Sequence:	GCTCTTCGGTTTCGAGGGCAAA
ACCN:	NM_006186 , NM_006186.1 , NM_006186.2 , NM_006186.3 , BC009288 , BC009288.2 , BC066890 , BM922519 , NM_006186.4
UniProt ID:	P43354
Synonyms:	HZF-3; NOT; NURR1; RNR1; TINUR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)