

Product datasheet for **HP209181**

PCBP1 Human qPCR Primer Pair (NM_006196)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5093
Forward Sequence:	GGACAACACACCATTTCTCCGC
Reverse Sequence:	AGCCTTTCACCTCTGGAGAGCT
ACCN:	<u>NM_006196</u> , <u>NM_006196.1</u> , <u>NM_006196.2</u> , <u>NM_006196.3</u> , <u>BC039742</u> , <u>BC039742.1</u> , <u>BM794597</u> , <u>NM_006196.4</u>
UniProt ID:	<u>Q15365</u>
Synonyms:	HEL-S-85; hnRNP-E1; hnRNP-X; HNRPE1; HNRPX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)