

Product datasheet for **HP209193**

PHKG1 Human qPCR Primer Pair (NM_006213)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5260
Forward Sequence:	CACAGACTTTGGCTTTTCCTGCC
Reverse Sequence:	ATGTCCACCTCTTCCCCTAGC
ACCN:	NM_006213 , NM_006213.1 , NM_006213.2 , NM_006213.3 , NM_006213.4 , BC051327 , BC051327.1 , BC069655 , BC069679 , BC069738 , BC069754 , BC074753 , NM_006213.5
UniProt ID:	Q16816
Synonyms:	PHKG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)