

Product datasheet for **HP209207**

POLE Human qPCR Primer Pair (NM_006231)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5426
Forward Sequence:	ACGCTGGAAGAGGTGTATGGCT
Reverse Sequence:	GGAACGGTTCTCAGAGATGAGC
ACCN:	<u>BC021559</u> , <u>NM_006231</u> , <u>NM_006231.1</u> , <u>NM_006231.2</u> , <u>NM_006231.3</u> , <u>BC021559.1</u> , <u>BC007599</u> , <u>BC011376</u> , <u>BC087613</u> , <u>BC144561</u> , <u>BC172436</u> , <u>BX647647</u> , <u>NM_006231.4</u>
UniProt ID:	<u>Q07864</u>
Synonyms:	CRCS12; FILS; IMAGEI; POLE1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)