

Product datasheet for **HP209213**

BRN3A (POU4F1) Human qPCR Primer Pair (NM_006237)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5457
Forward Sequence:	AGTACCCGTCGCTGCACTCCA
Reverse Sequence:	TTGCCCTGGGACACGGCGATG
ACCN:	NM_006237 , NM_006237.1 , NM_006237.2 , NM_006237.3 , BC148330 , BC148792 , BM714015 , BM726529 , NM_006237.4
UniProt ID:	Q01851
Synonyms:	ATITHS; brn-3A; BRN3A; Oct-T1; RDC-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)