

Product datasheet for **HP209230**

PKC delta (PRKCD) Human qPCR Primer Pair (NM_006254)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5580
Forward Sequence:	GCTGACACTTGCCGCAGAGAAT
Reverse Sequence:	GCCTTTGCCTGGATGTGGTAC
ACCN:	<u>NM_006254</u> , <u>NM_006254.1</u> , <u>NM_006254.2</u> , <u>NM_006254.3</u> , <u>BC043350</u> , <u>BC043350.1</u> , <u>BC127261</u> , <u>BE048101</u> , <u>NM_006254.4</u>
UniProt ID:	<u>Q05655</u>
Synonyms:	ALPS3; CVID9; MAY1; nPKC-delta; PKCD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)