

Product datasheet for **HP209249**

SRSF6 Human qPCR Primer Pair (NM_006275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6431
Forward Sequence:	GTGCTTTGGACAACTGGATGGC
Reverse Sequence:	CTCCTACTTCGTGACCGTCTTC
ACCN:	NM_006275 , NM_006275.1 , NM_006275.2 , NM_006275.3 , NM_006275.4 , NM_006275.5 , BC006832 , BC006832.2 , BM788391 , BM830876
UniProt ID:	Q13247
Synonyms:	B52; HEL-S-91; SFRS6; SRP55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)