

Product datasheet for **HP209259**

CD90 (THY1) Human qPCR Primer Pair (NM_006288)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7070
Forward Sequence:	GAAGGTCCTCTACTTATCCGCC
Reverse Sequence:	TGATGCCCTCACACTTGACCAG
ACCN:	<u>NM_006288</u> , <u>NM_006288.1</u> , <u>NM_006288.2</u> , <u>NM_006288.3</u> , <u>BC065559</u> , <u>BC005175</u> , <u>BC037832</u> , <u>BE672723</u> , <u>BM504140</u> , <u>NM_006288.5</u>
UniProt ID:	<u>P04216</u>
Synonyms:	CD90; CDw90
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)