

## Product datasheet for **HP209261**

### TNFAIP3 Human qPCR Primer Pair (NM\_006290)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7128
Forward Sequence:	CTCAACTGGTGTGCGAGAAGTCC
Reverse Sequence:	TTCCTTGAGCGTGCTGAACAGC
ACCN:	<a href="#">NM_006290</a> , <a href="#">NM_006290.1</a> , <a href="#">NM_006290.2</a> , <a href="#">NM_006290.3</a> , <a href="#">BC041790</a> , <a href="#">BC064689</a> , <a href="#">BC098379</a> , <a href="#">BC113871</a> , <a href="#">BC114480</a> , <a href="#">BQ002744</a> , <a href="#">NM_006290.4</a>
UniProt ID:	<a href="#">P21580</a>
Synonyms:	A20; AISBL; OTUD7C; TNFA1P2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)