

Product datasheet for **HP209265**

UQCRB Human qPCR Primer Pair (NM_006294)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7381
Forward Sequence:	TAAGAGGGCACTGGACCTGAAC
Reverse Sequence:	CATGATTACTTCTTTGCCATTCTT
ACCN:	NM_006294 , NM_006294.1 , NM_006294.2 , NM_006294.3 , NM_006294.4 , BC005230 , BC005230.1 , BC041005 , BI547559 , BU535281 , NM_006294.5
UniProt ID:	P14927
Synonyms:	MC3DN3; QCR7; QP-C; QPC; UQBC; UQBP; UQCR6; UQPC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)