

Product datasheet for **HP209275**

SEM1 Human qPCR Primer Pair (NM_006304)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7979
Forward Sequence:	CTGGGCTGGCTTAGATGAAGATG
Reverse Sequence:	CTTCTGGATGCTATGAAGTCTCC
ACCN:	NM_006304 , NM_006304.1 , BC032782 , BC032782.1
UniProt ID:	P60896
Synonyms:	C7orf76; DSS1; ECD; PSMD15; SHFD1; Shfdg1; SHFM1; SHSF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)