

Product datasheet for **HP209283**

USP15 Human qPCR Primer Pair (NM_006313)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9958
Forward Sequence:	CCCTACACTGCTGTAAGGACCA
Reverse Sequence:	GAGGGAAGTTCTTGATCCTGGC
ACCN:	<u>NM_006313, NM_006313.1, BC125123, BC020688, BC063454, BC110908, NM_006313.3</u>
UniProt ID:	<u>Q9Y4E8</u>
Synonyms:	UNPH-2; UNPH4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)