

Product datasheet for **HP209287**

BASP1 Human qPCR Primer Pair (NM_006317)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10409
Forward Sequence:	AAGCGCCTAGTTCCACACCCAA
Reverse Sequence:	CTGTCCTTGCTCACTCTTTCACGG
ACCN:	NM_006317 , NM_006317.1 , NM_006317.2 , NM_006317.3 , NM_006317.4 , BC000518 , BG704783 , BM888431 , BM973114 , BT019340 , BT019341 , BU181147 , NM_006317.5
UniProt ID:	P80723
Synonyms:	CAP-23; CAP23; NAP-22; NAP22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)