

Product datasheet for **HP209327**

SEC23 (SEC23A) Human qPCR Primer Pair (NM_006364)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10484
Forward Sequence:	TACAGCAGCCACCTCCTTCCAA
Reverse Sequence:	CGCAAAGGTCTCTTCCCTGTG
ACCN:	NM_006364 , NM_006364.1 , NM_006364.2 , NM_006364.3 , BC036649 , BC036649.1 , BM541968 , BX537566 , NM_006364.4
UniProt ID:	Q15436
Synonyms:	CLSD; hSec23A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)