

Product datasheet for **HP209333**

CRTAP Human qPCR Primer Pair (NM_006371)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10491
Forward Sequence:	GTCATGCAGCAGAACCTGGTGT
Reverse Sequence:	CTCCTTCTGGAGTGTGGTCACA
ACCN:	NM_006371 , NM_006371.1 , NM_006371.2 , NM_006371.3 , NM_006371.4 , BC008745 , BC008745.2 , BC073904 , BC094866 , BX648233 , NM_006371.5
UniProt ID:	O75718
Synonyms:	CASP; LEPREL3; OI7; P3H5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)