

Product datasheet for **HP209339**

Semaphorin 4D (SEMA4D) Human qPCR Primer Pair (NM_006378)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10507
Forward Sequence:	AGCTCTGCACAAAGCCATCAGC
Reverse Sequence:	CCAGCATAGACAAACCTGTTGCC
ACCN:	<u>NM_006378</u> , <u>NM_006378.1</u> , <u>NM_006378.2</u> , <u>NM_006378.3</u> , <u>BC054500</u> , <u>BC054500.1</u> , <u>BC137515</u> , <u>BC137518</u> , <u>BX648216</u>
UniProt ID:	<u>Q92854</u>
Synonyms:	A8; BB18; C9orf164; CD100; coll-4; COLL4; GR3; M-sema-G; SEMAJ
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)