

## Product datasheet for **HP209356**

### Ribonuclease H2, subunit A (RNASEH2A) Human qPCR Primer Pair (NM\_006397)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10535
Forward Sequence:	GCCGTGAAGAAATGGCAGTTCG
Reverse Sequence:	GTGCTCCTCAACCACGCTTTTG
ACCN:	<a href="#">NM_006397</a> , <a href="#">NM_006397.1</a> , <a href="#">NM_006397.2</a> , <a href="#">BC011748</a> , <a href="#">BE560339</a> , <a href="#">BG825340</a> , <a href="#">BM664190</a> , <a href="#">NM_006397.3</a>
UniProt ID:	<a href="#">O75792</a>
Synonyms:	AGS4; JUNB; RNASEHI; RNHIA; RNHL; THSD8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)