

## Product datasheet for **HP209358**

### BATF Human qPCR Primer Pair (NM\_006399)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10538
Forward Sequence:	GATGTGAGAAGAGTTCAGAGGAG
Reverse Sequence:	GTTTCTCCAGGTCTTCGCTCTC
ACCN:	<u><a href="#">NM_006399</a></u> , <u><a href="#">NM_006399.1</a></u> , <u><a href="#">NM_006399.2</a></u> , <u><a href="#">NM_006399.3</a></u> , <u><a href="#">BC032294</a></u> , <u><a href="#">BC032294.1</a></u> , <u><a href="#">BP235761</a></u> , <u><a href="#">NM_006399.5</a></u>
UniProt ID:	<u><a href="#">Q16520</a></u>
Synonyms:	B-ATF; BATF1; SFA-2; SFA2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)