

## Product datasheet for **HP209367**

### ARPC1A Human qPCR Primer Pair (NM\_006409)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10552
Forward Sequence:	AGAGTTCCTGCCGCTCCTAAGT
Reverse Sequence:	ACTTGAGACGAAGGTCAGGCA
ACCN:	<u><a href="#">NM_006409</a></u> , <u><a href="#">NM_006409.1</a></u> , <u><a href="#">NM_006409.2</a></u> , <u><a href="#">NM_006409.3</a></u> , <u><a href="#">BC039594</a></u> , <u><a href="#">BC039594.1</a></u> , <u><a href="#">BC006445</a></u> , <u><a href="#">BC047889</a></u> , <u><a href="#">BC054027</a></u> , <u><a href="#">BM924483</a></u> , <u><a href="#">NM_006409.4</a></u>
UniProt ID:	<u><a href="#">Q92747</a></u>
Synonyms:	Arc40; HEL-68; HEL-S-307; SOP2Hs; SOP2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)