

## Product datasheet for **HP209383**

### TCP1 delta (CCT4) Human qPCR Primer Pair (NM\_006430)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10575
Forward Sequence:	CATGTCTCGACCTGTGGAAGCTG
Reverse Sequence:	CACTCATTGGAGAAAGCAGACTTG
ACCN:	<a href="#">NM_006430</a> , <a href="#">NM_006430.1</a> , <a href="#">NM_006430.2</a> , <a href="#">NM_006430.3</a> , <a href="#">BC014676</a> , <a href="#">BC093446</a> , <a href="#">BC106933</a> , <a href="#">BC106934</a> , <a href="#">BM765256</a> , <a href="#">NM_006430.4</a>
UniProt ID:	<a href="#">P50991</a>
Synonyms:	CCT-DELTA; Cctd; SRB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)