

Product datasheet for **HP209391**

TXNRD2 Human qPCR Primer Pair (NM_006440)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10587
Forward Sequence:	GCACCTTTGACACCGTCCTGTG
Reverse Sequence:	CACCAGGATCTTCTGAGTGTCG
ACCN:	<u>NM_006440</u> , <u>NM_006440.1</u> , <u>NM_006440.2</u> , <u>NM_006440.3</u> , <u>BC007489</u> , <u>BC032029</u> , <u>BC038708</u> , <u>BC099923</u> , <u>BC117354</u> , <u>BX957216</u> , <u>NM_006440.5</u>
UniProt ID:	<u>Q9NNW7</u>
Synonyms:	GCCD5; SELZ; TR; TR-BETA; TR3; TRXR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)