

Product datasheet for **HP209392**

MTHFS Human qPCR Primer Pair (NM_006441)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10588
Forward Sequence:	TGGAATATCCCTCAGCCTGGTG
Reverse Sequence:	CCAGTCGGTTGCCATGTTTGTC
ACCN:	NM_006441 , NM_006441.1 , NM_006441.2 , NM_006441.3 , BC019921 , BC019921.1 , BC012417 , BC037852 , BM313547 , BM450865 , NM_006441.4
UniProt ID:	P49914
Synonyms:	HsT19268; NEDMEHM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)