

Product datasheet for **HP209396**

OATP2 (SLCO1B1) Human qPCR Primer Pair (NM_006446)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10599
Forward Sequence:	CGTAGAGCAACAGTATGGTCAGC
Reverse Sequence:	TTGGCAATTCCAACGGGTGTCAG
ACCN:	<u>NM_006446</u> , <u>NM_006446.1</u> , <u>NM_006446.2</u> , <u>NM_006446.3</u> , <u>NM_006446.4</u> , <u>BC005243</u> , <u>BC020696</u> , <u>BC054878</u> , <u>BC070172</u> , <u>BC105596</u> , <u>BC114376</u>
UniProt ID:	<u>Q9Y6L6</u>
Synonyms:	HBLRR; LST-1; LST1; OATP-C; OATP1B1; OATP2; OATPC; SLC21A6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)