

## Product datasheet for **HP209433**

### Tbp7 (PSMC4) Human qPCR Primer Pair (NM\_006503)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5704
Forward Sequence:	GCAAGATGAACCTCTCTGAGGAG
Reverse Sequence:	CGGTTTTACGGACAGCCAACA
ACCN:	<u><a href="#">NM_006503</a></u> , <u><a href="#">NM_006503.1</a></u> , <u><a href="#">NM_006503.2</a></u> , <u><a href="#">NM_006503.3</a></u> , <u><a href="#">BC010396</a></u> , <u><a href="#">BC014488</a></u> , <u><a href="#">BC014488.1</a></u> , <u><a href="#">BC000343</a></u> , <u><a href="#">BC018811</a></u> , <u><a href="#">BI828251</a></u> , <u><a href="#">BT007232</a></u> , <u><a href="#">NM_006503.4</a></u>
UniProt ID:	<u><a href="#">P43686</a></u>
Synonyms:	MIP224; RPT3; S6; TBP-7; TBP7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)