

Product datasheet for **HP209443**

Seryl tRNA synthetase (SARS) Human qPCR Primer Pair (NM_006513)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6301
Forward Sequence:	CCCAGAGAATGTGCTGAGTTTCG
Reverse Sequence:	CTCAAACCGCTCTGCTTCCAAC
ACCN:	<u>NM_006513</u> , <u>NM_006513.2</u> , <u>NM_006513.3</u> , <u>BC009390</u> , <u>BC009390.2</u> , <u>BC000716</u> , <u>BC111398</u> , <u>BM666450</u> , <u>BM798814</u> , <u>BP363347</u> , <u>NM_006513.4</u>
UniProt ID:	<u>P49591</u>
Synonyms:	NEDMAS; SARS; SERRS; SERS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)